


ORDER OF LAB VALUE DRAWS



ADDITIVES BY TUBE COLOR

Sterile - SPS

Lt Blue - Sodium Citrate

Plain Red - No Additive

Plastic Red - Clot Activator

Gold, Red/Gray - SST/Gel w/Clot Activator

Green - Heparin-Lithium or Sodium

Lavender, Pink - EDTA

Gray - Sodium Fluoride, Potassium Oxalate

OF TIMES TO INVERT TUBE

Clot Activator - 5 times

Sodium Citrate - 3-4 times

Heparin - 8 times

EDTA - 8 times

Sodium Fluoride - 8 times